

Systeme für Aquakultur,  
Aquaristik, Labore und  
zur Wasseraufbereitung

Systems for aqua culture,  
sea water aquaria, labs and  
water desalination and purification


Systèmes pour aquaculture,  
aquariums eau de mer,  
laboratoires et traitements d'eau




AquaCare GmbH & Co. KG  
Am Wiesenbusch 11  
D-45966 Gladbeck  
Tel.: +49-2043-375758-0  
Fax: +49-2043-375758-90  
<http://www.aquacare.de>  
e-mail: info@aquacare.de

# KWR70 ... KWR300

## KWR70 & KWR100




# KWR150 ... KWR300



KWR150: FK-M08x25-PA  
 KWR200: FK-M08x25-PA  
 KWR300: SS-M10x030-PA



KWR150: U8,5-PA  
 KWR200: U8,5 PA  
 KWR300: U11 PA



PA-W1/4-4

KWR150: PA-W1/2-12

PMMA-cap-R-150  
 PMMA-cap-R-200  
 PMMA-cap-R-300


KWR150: 160x4N  
 KWR200:  
 KWR300:

KWR150: 004N70  
 KWR200:  
 KWR300:

Si07-04

Si24-20

RE2368



sieb-150-KWR  
 sieb-200-KWR  
 sieb-300-KWR




KWR150: Pumpenrad1048  
 Eheim #: 7645990  
 KWR200+300: Pumpenrad1250  
 Eheim #: 7640900



sieb-150-4mm  
 sieb-200-4mm  
 sieb-300-4mm

SK-M08x25-PA

SK-M06x16-PA



fuss22x11

PPSV010822W


PA 1/4"-4

KWR150: Eheim 1048  
 KWR200: Eheim 1250  
 KWR300: Eheim 1250








KWR150: 730-025  
 KWR200: 730-025  
 KWR300: 730-032


SK-M08x25-PA

KWR150: KWR-H150  
 KWR200: KWR-H200  
 KWR300: KWR-H300




RO/VE  
RO/DI





6.



7.

